# dima® Print Denture Base Try-in

dima® Print Denture Base Try-in

- DE) Gebrauchsanweisung
- GB) Instructions for use
- FR) Mode d'emploi
- ES) Instrucciones de uso
- IT) Istruzioni per l'uso
- PT) Instruções de uso
- NL) Gebruiksaanwijzing

**Explanation of symbols on labelling**

- MD Medical device
- LOT Batch code
- Manufacturer
- Importer
- Keep away from sunlight
- Consult instructions for use
- Use-by date
- Storage temperature

Manufacturer:  
DENTCA, INC.  
357 VAN NESS WAY, STE 250  
TORRANCE, CA 90501  
USA

EU Importer:  
Kulzer GmbH  
Leipziger Straße 2  
63450 Hanau (Germany)  
+49 800 4372522

Distributed in USA /  
Canada exclusively by:  
Kulzer, LLC  
4315 South Lafayette Blvd.  
South Bend, IN 46614-2517  
1-800-431-1785

Caution: Federal law restricts  
this device to sale by or on the  
order of a dental professional.

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EU Representative:  
MT Promed Consulting GmbH  
Ernst-Heckel-Strabe 7  
66386 St. Ingbert  
Germany

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## GB) Instructions for use dima® Print Denture Base Try-in

**Indications for Use**  
This product is light-curable resin indicated for fabrication and repair of removable dentures and baseplates. The material is an alternative to traditional heat-curable and auto polymerizing resins.  
Fabrication of dental prosthetics with this product requires a computer-aided design and manufacturing (CAD/CAM) system that includes the following components: digital denture files based on a digital impression, stereolithographic additive printer, and curing light equipment.  
This product is particularly suitable to fabricate the try-in denture before installing a final denture.

- Requirements**
1. Digital denture model file with minimum thickness 2 mm; STL format file
  2. Additive printer and its operation software;

Additive Printer	Operation Software	Provider
cara Print 4.0 or cara Print 4.0 pro	cara Print CAM or cara Print CAM 2.0	Kulzer
Asiga Max, Pro2 or Pro 4K	Asiga Composer	Asiga

3. Curing light equipment

Curing Equipment	Provider
HiLite power 3D / LEDcure	Kulzer

### Specific Manufacturing Considerations

1. Digital denture model file
  - 1.1 File format: STL file
  - 1.2 Minimum thickness of denture model: >= 2 mm
  - 1.3 File size should be uploadable in the 3D printer operation software.
2. Additive printer
  - 2.1 Hardware
    - a. Wavelength : 385 nm or 405 nm
    - b. Light source
      - Stereolithographic (SLA) method; laser with 25 mW < X < 250 mW
      - Digital Light Processing (DLP) method; high power LED or lasers
    - c. Build Volume: > 103 x 58 x 130 mm (Least fit one arch)
    - d. Build Path: line drawing path or surface layer drawing path
  - 2.2 Features of Operation Software
    - a. STL file import
    - b. Automatic rotation and placement
    - c. Layer slicer for path inspection
    - d. Auto and manual generation of supports
  - 2.3 Printing Parameter

Printer Model	Layer Thickness (micron)	Recommended orientation angle (degree)	Support point size (mm)	Support density (mm)
cara Print 4.0 or cara Print 4.0 pro	50, 100	40–60	1–2	min. 1
Asiga Max, Pro2 or Pro 4K	50–100	20–90	1.2–1.5	spacing 5.0 mm

- 2.4 Environmental Conditions
  - a. Temperature: 18–30°C (64–86°F)
  - b. Relative Humidity: 30–90%
- 2.5 Cleaning Kit
  - a. Rinse bath and tubs, flush cutter, paper towel, squeeze bottle for isopropyl alcohol, Scraper

3. Recommended Curing light equipment (Post curing units)

3.1 Flood Type Curing Equipment

Manufacturer / Model	Curing Chamber	Supply voltage (V / Hz)	Lamp power (W)	Light intensity (mW/cm²)	Lamp wavelength (nm)	Curing time (min)
Kulzer / HiLite power 3D	Required	100, 115, 230 / 50–60	200	N/A	390–540	2 x 10
Kulzer / LEDcure	Required	100–240 / 50, 60	150	N/A	370–470	select dima Try-in program

- 3.2 Accessories
- a. USP Grade glycerin
  - b. Transparent glass container and 2 glass plates

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Stand: 2023-06

- c. Heat-protective gloves and silicone coated stainless steel tong
  - d. Thermocouple
4. Notification
- 4.1 The device specifications have been validated using the software, printers, and process parameters specified in this document. Any other printers, operation software and post-printing processes will be outside of the device specifications and the FDA clearance. Users shall follow this document to use the device.
  - 4.2 If there is a serious incident (death or permanent damage to a patient) that has occurred in relation with this device please report to DENTCA (info@dentca.com) or your local authority of medical devices.
  - 4.3 When you receive damaged or unintentionally open bottles before use, or if the packaging is exposed to environmental conditions outside of the specified in the label, please inform to DENTCA (info@dentca.com).

### Warnings:

1. dima Print Denture Base Try-in resin contains polymerizable monomers which may cause skin irritation (allergic contact dermatitis) or other allergic reactions in susceptible persons. If contact with skin, wash thoroughly with soap and water. If skin sensitization occurs, discontinue use. If dermatitis or other symptoms persist, seek medical assistance.
2. Avoid inhalation or ingestion. High concentration vapors may cause headache, irritation of eyes or respiratory system. Direct contact with eyes may cause possible corneal damage. Long-term excessive exposure to the material may cause more serious health effects. Monitor air quality per OSHA standards.
3. Eye Contact: Immediately flush eyes with plenty of clean water for at least 20 minutes, and consult a physician. Wash the contacted area thoroughly with soap and water. Inhalation: In case of exposure to a high concentration of vapor or mist, remove person to fresh air. Give oxygen or artificial respiration as required.
4. Ingestion: Contact your regional poison control center immediately.
5. BURN HAZARD: GLYCEROL BATH CAN REACH TEMPERATURES OF 90°C (-200°F) AND LEAD TO SEVERE BURNS. Only trained users should perform the glycerol curing step with caution and appropriate PPE. We also recommend placing a warning label on the window of the cure unit to alert all lab users to the potential hazard.

### Precautions:

1. When washing the printed Try-in denture with isopropanol or grinding the Try-in denture, it should be in a properly ventilated environment with proper protective masks and gloves.
  2. Store dima Print Denture Base Try-in resin at 15–25°C (60–77°F) and avoid direct sunlight. Keep container closed when it is not in use. Product shall not be used after expiration date.
  3. Expired or unused dima Print Denture Base Try-in should be completely cured or polymerized prior disposal.
- Adverse Reactions:**
1. Direct contact with the uncured resin may induce skin sensitization in susceptible individuals.
  2. Proper ventilation and personal protective equipment should be used when grinding printed Try-in denture as the particulate generated during grinding may cause respiratory, skin and eye irritation.

### Procedure to Fabricate the Try-in Denture

1. Printing Preparation
  - a. Open the 3D printer cover and fill the resin tank or tray of the printer with dima Print Denture Base Try-in resin up to the required filling line by manufacturer. (When filling the resin into the resin tank or tray, gloves and mask should be used.)
  - b. Close the printer cover.
2. Printing
  - a. Load the denture model file to be printed in printer operation software which printer manufacturer recommended.
  - b. Use auto-orientation or manual orientation to find its optimal position for printing. The recommended orientation by printer provider is a tilted orientation such as space diagonal from 40–60° angle. If auto-orientation is not satisfied, rotate to make optimal position.
  - c. Generate support sticks on the denture model. The support is not enough, add supports on the denture model.
  - d. Use layout tools of the software to move the denture model within the imaginary build platform to prevent the overlapping between models.
  - e. Start printing.
3. Cleaning
  - a. Detach the printed Try-in denture from the build platform.
  - b. Use a small flush cutter to remove the support sticks from the printed Try-in denture.
  - c. Wash the Try-in denture with isopropyl alcohol.
  - d. Use air blowing to dry the Try-in denture or dry it at room temperature under ventilation system or open area.
4. Post Processing of Try-in denture
  - a. Add dima Print Denture Base Try-in material manually on the lingual ridge area to thicken and cure it. It is recommended to make it thicker than 4.0mm for maxillary and 4.5mm for mandibular.
  - b. Cure the printed Try-in denture by sinking into glycerin container for the required curing time under recommended post-curing unit. **Note:** It is recommended the glycerol to replace every 80 hrs running or every three months whichever comes first.
  - c. Take out the printed denture from the curing oven using coated tong (**Be careful hot glycerin**).
  - d. Rinse the cured denture with a water.
  - e. Smooth the support marks using a bur after washing the Try-in denture with water and drying.
  - f. Please make sure all the support marks completely grind.

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