
Xthetic®
denture Fx

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(GB) Instructions for use

Product Description
Resin for customization of full and partial dentures. Compatible with all common PMMA-based denture materials.

Composition

Two-component acrylic resin consisting of powder and liquid. Powder main component: Methyl methacrylate copolymer. Liquid main components: Methyl methacrylate, dimethacrylate.

Processing in the casting technique

The colors are applied direct with a natural-hair brush to the desired position of the duplicating gel or silicone matrix.

- Customization of full and partial dentures using the cuvette or matrix technique
- Subsequent customization of finished full and partial dentures

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Notes

- The intensive colors can be mixed with one another as well as with all common PMMA-based denture materials.

Subsequent customization of the denture

- For processing the material it is recommended to use a tray on which the colors can be spread and mixed.

Composition

- When using the material, the modeling liquid supplied should be used to achieve the longest possible working time and to prevent premature drying of the material.

Polymerization

Polymerization takes place under a pressure of 2-3 bar and a water temperature of 45-50°C for approx. 15 minutes.

After completion

The denture must be stored in water until incorporation. It is generally recommended to store the denture in water for 24 hours to further reduce the residual monomer content (< 4.5% according to ISO 20795-1) and thus prevent the risk of irritation.

Shades

The material is available in gingiva light, gingiva dark, ivory, rose, indian, crimson, orchid, sienna, maroon and indigo

Forms of delivery

Powder: 5g, 35g, 100g
Liquid: 80ml
Set pack: contains all colors, modeling liquid and mixing cup.

® = registered trademark of AcrylX® GmbH

Storage and Safety Measures

Do not use the material after the expiration date. Store at temperatures not exceeding 25°C [77°F]. Protect from sunlight. Close containers

immediately after use. Keep away from heat/sparks/open flames/hot surfaces. No smoking. Dispose of contents/container in accordance with local disposal regulations.

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® = marca registrada de AcrylX® GmbH
® = marchio registrato di AcrylX® GmbH