

JULY / AUGUST 2025 HIGHLIGHTS

durodent dental

DIGITAL CAD / CAM PRODUCT RANGE



ceraMotion Z "Hybrid" Zirconia 14mm / 18mm / 22mm

Multi-shaded and multi-layered zirconium oxide made of two different zirconium oxide materials for all restorations in the anterior and posterior area.
Suitable for all open CAD/CAM systems - 98mm Dia.

14MM FROM

\$238.57 ea

ex GST

INTRO OFFER ***
10% OFF

*** Some ranges and sizes are SPECIAL ORDER 2-3 week leadtime



PMMA Milling Discs For Premium Denture Base and Teeth



Shaded PMMA (MONO) Tooth Shades

Full Shade Range

Achieve great results
at a Value for Money
Price.



20mm x 98mm

\$67.85

ex GST

BUY 3 + SAVE 10%

Multi-Layer PMMA Tooth Shades

Full Shade Range
13 Intricate Layers

Resulting in **highly
esthetic outcomes**, dentin
and enamel blended
seamlessly.



20mm x 98mm

\$111.14

ex GST

BUY 3 + SAVE 10%

Lucitone Denture Base High Impact PMMA

30mm or 25mm x
98mm Dia

Original Pink or
Light Pink



30mm x 98mm

\$104.52

ex GST

BUY 3 + SAVE 10%

SCAN SPRAYS



Scantist 3D Vanishing
(Self Evaporating)



Ernst Hinrich
Scan Spray PLUS

SOFT SPLINT CAD DISC



Reliance
Dental USA
DuraSprint
98mm x
20mm

20mm x 98mm

\$125.00

ex GST

BUY 3 + SAVE 10%

3D PRINT RESINS



SEE RESIN
RANGES
ONLINE



Dentsply
Sirona

Subscribe @ Durodent.com.au for our monthly email broadcast
Keep informed – New products, highlights, special offers

SUBSCRIBE

Valid JULY / AUGUST 2025 - All Prices Exclude GST - Prices Do Not Include Freight - Durodent Reserves The Right To Amend Errors and Prices - Subject to stock availability

ONLINE STORE
www.durodent.com.au

SAVE TIME - ORDER ONLINE



Speak to our team
03 9720 6700

LAB PUTTY

Detax Blue Eco Lab Putty 1:1 Mix

Equal Part Mix 1:1
Great for processing dentures
Final Hardness: 86 Shore A.
Heat Resistant up to 200 °C.



1.6kg (2 x 800g)

\$89.10 ex GST



10kg (2 x 5kg)

\$369.00 ex GST

10% OFF

IMPRESSION DISINFECTANT

Cavex ImpreSafe Starter Kit HA053

3 Minute soak protocol -
To kill bacteria, fungi and
viruses.

33 Weeks From 1L Bottle

Includes: 1 x 1 litre Cavex
ImpreSafe Disinfection Fluid, 1 x
Disinfection Container, 1 x Timer,
50 x Plastic Bags



STARTER KIT

\$218.25 ex GST

10% OFF

MODEL TRIMMER

Renfert MT3 with Marathon Diamond Disc (Partial) 18080500

WET TRIMMER

Top performing industry
favourite. Easy to clean and
maintain.

*Renfert 100 Year Anniversary
Extended 3+1 Year Guarantee
Offer - Register with Renfert*



\$2450.00 ex GST



DENTURE POLISH

Ainsworth Dentishine

Dental polishing paste.

Used to obtain a high shine
on a finished acrylic denture
with soft buff.

250G SMALL

\$20.25 ex GST

500G LARGE

\$29.25 ex GST



POLISHING

Renfert Ceramic Restorations ALL IN ONE POLISH Introduction Starter Kit 5100103

18g All In One
Wide spectrum of
diamond grain sizes
+ Instruction Guide
+ HP Brushes and Buff



SPECIAL

\$65.00 ex GST



BITE REGISTRATION

Zhermack Occlufast + Color 2 x 50ml C200781

Thermochromic A-
Silicone

Clinical Use: Bite
Registration

Laboratory Uses –
Insulating / supporting
around teeth during
processing / injection
technique.



BUY 3 SAVE 10%

\$85.00 ea ex GST

BUY 6 GET 1 FREE

Valid JULY / AUGUST 2025 - All Prices Exclude GST - Prices Do Not Include Freight - Durodent Reserves The Right To Amend Errors and Prices - Subject to stock availability

JULY / AUGUST 2025 HIGHLIGHTS

durodent
dental

Kulzer Venus White Pro 16% Carbamide Peroxide Tooth Whitening

10% OFF

Mint flavoured gel whitening system with sensitivity control.

Extended wear or overnight application.



PATIENT KIT 6 x 1.2ml

\$65.41

ex GST

REFILL PACK 3 x 1.2ml

\$37.02

ex GST

Ultradent LC Block Out Resin 1.2ml Syringe

10% OFF

A low viscosity, blue pigmented light-cured resin.

Provides reservoir space for bleaching trays.



4 PACK

\$88.20

ex GST

BULK PACK 20PC

\$320.40

ex GST

"Zeta Med" Snoring Device Component With Wings – Type C - Titra Table Appliance - 1 Unit

10% OFF

Snoring Device Appliance component for lab use manufacturing a custom patient snoring device.



\$120.60

ex GST

Cello Econaps Bibs - White 4 Ply - 1000/Ctn

10% OFF

White Bib Recycled Paper With Plastic Barrier

208 x 280mm (AEN)



\$76.50

ex GST

AcrylX Xthetic Temp Tooth Shade Acrylic Starter Kit (1-608-01) - Made In Germany

10% OFF

High Grade Strong
C&B Denture Tooth
Colour Resin

Colour Stable

Includes: 6 x 35g
Powders (A1, A2, A3,
A3.5, A4, B3)
4 x 5g (Clear, IC1, IC21
IC3), 1 x 100ml Liquid,
1 x Mixing Cup



\$337.50

ex GST

Kulzer Palabond Tooth to Base Bond - 45ml - 64708082

10% OFF

Increases the bond between
**denture teeth and denture
base**, as well as the bond
between repair material and
base material.

**** 10 Minute Maximum
working time. Complete
process within 10 minutes**



\$54.09

ex GST



ADD TO CART

ONLINE ORDER FREIGHT OFFERS

Free Shipping or Subsidised Flat Rate Shipping Offers Available
for most orders valued over \$250.00 ex GST
with 2 or more items when placed Online through our Web Shop

*Conditions Apply

Excludes some items (eg. Stone, Equipment, Large Items) and delivery locations.
It will tell you at checkout if promo applies.
Freight promo not available for phone orders.

ONLINE STORE
www.durodent.com.au

SAVE TIME - ORDER ONLINE

ADD TO CART



Speak to our team
03 9720 6700



Durodent welcomes the Dentalife Infection Control range
NOW IN STOCK - WIPES - SPRAY - MOUTHRINSE



**Clinicare Neutral
Detergent Range**



**Clinicare "Alcohol Free"
Disinfectant Range**



**Clinicare HGD Hospital
Grade Disinfectant Range**



**Swirl Mouthrinse
Alcohol Free**



**Swirl PREP
"Pre-Procedural"
Disinfectant
Mouthrinse**

10% OFF DENTALIFE RANGE - INTRODUCTION OFFER



Liquid Form



Canister Refill Packs



Dispensers



Canister Wipes



Flat Pack Wipes

**FROM ONLY
\$13.50 Per Pack ex GST**

LAB YELLOWSTONES AND DIESTONES

Durodent PROS STONE YELLOWSTONE

The Perfect Lab Stone

PSI 8,500 (58Mpa)
Set Time: 11-14 Mins
Exp: 0.12%
Ratio: 30ml/100g
Made In USA



5kg Pail
\$22.50 ex GST

20kg Box
\$80.10 ex GST

Durodent CHROME STONE DIESTONE - Blue

Great Chrome Model Stone

PSI 13,000 (90Mpa) -
Set Time: 12 - 15 Mins -
Exp: 0.15% - Ratio: 24ml/100g
Made In USA



5kg Pail
\$40.50 ex GST

20kg Box
\$108.00 ex GST

Ernst Hinrichs LABORIT Yellowstone Buff 20kg Box

Popular all round yellow stone

PSI 8,500 (58Mpa) - Set Time: 10-12 Mins - Exp: 0.14% -
Ratio: 30ml/100g
Made In Germany



20kg Box
\$71.91 ex GST

Kulzer Hera MOLDASTONE CN, Light Brown, 18kg (4 x 4.5kg) 66014745

*High Strength Dental Stone,
Type 4. Hard dental diestones for
fixed prosthesis.*

PSI 9,500 (65Mpa) - Set Time: 12 Mins - Exp: 0.08% - Ratio: 22ml/100g



18kg Box
\$130.50 ex GST

Are you interested in trying the
proven top selling Kulzer Pala Denture Teeth range???



Incredible VALUE Introduction Kits On Offer

Low risk, great opportunity to give the teeth a really good go
and experience the difference.

Kulzer DELARA 6/8
Small Assortment Intro Case
66077520

**PROVIDES A MASSIVE 38%
DISCOUNT OFF RRP**

RESTRICTED OFFER TERMS *

\$780.00 ex GST

Includes Shades A2/A3
48 Anteriors / 36 Posteriors
Shade Guide
All Moulds



**Kulzer PALA MONDIAL 6/8
and PREMIUM 6**

Small Assortment Intro Cases

48 Anteriors / 36 Posteriors - 2 Draws
Selection Of All Mould To Try
Two Shades - Options available

**PROVIDES A MASSIVE 55%
DISCOUNT OFF RRP***

MONDIAL 6/8
\$1484.00 ex GST

PREMIUM 6 / MONDIAL 8
\$1814.00 ex GST




*RESTRICTED OFFER SORRY 😊 - 1 KIT PER CUSTOMER / ACCOUNT ONLY - IT IS A VERY SPECIAL OFFER TO HELP YOU GET STARTED.

IF YOU HAVE ALREADY PURCHASED ONE OF THESE IN THE PAST, WE WILL HAVE TO KINDLY SAY NO, SORRY.

PLEASE HAVE A LOOK AT THE OTHER GREAT VALUE CUSTOM SETUP OFFERS AS THE NEXT LEVEL.

Want to use **CARDED TEETH** in your **DIGITAL WORKFLOW?**

*Digital Tooth Library licence subscription maybe required



	Premium	Mondial	Idealis	Delara
	Pala Mix & Match System			
3shape	✓	✓	✓	✓
exocad		✓	✓	✓

3shape

exocad

KULZER
MITSUI CHEMICALS GROUP

ONLINE STORE
www.durodent.com.au

SAVE TIME - ORDER ONLINE



Speak to our team
03 9720 6700

DURODENT ACRYLIC TEETH RANGE

durodent
dental



KULZER
MITSUI CHEMICALS GROUP

DELARA 3 Layer Modern Aesthetics

Everybody's Delara – a modern all-round tooth line for everyday use with lifelike aesthetics.

A single tooth line that meets every need of dental laboratories – delivering aesthetics, form and function, with a lively surface structure and unique layering, along with modern and bright shading for a fresh and vivid appearance.

3 Layer Incomp Cad Cam

Shades Available: A1 – D4 + 2 Bleach Shades BL2 / BL3

Moulds: 18 Upper / 5 Lower Anterior Moulds - 8 Posteriors Moulds



Only
\$14.85
Per Card or Less

BUY 50-99
SAVE 10%
\$13.36 ea

BUY 100+
SAVE 15%
\$12.62 ea

INTRO CASE +
CUSTOM BULK
OFFERS
AVAILABLE

Kulzer Pala MONDIAL Acrylic Teeth – CAD CAM Designed Perfection – Wear Resistance Tech

Superior Wear Resistance Properties For Fixed or Removable Prosthetics

Highly Pleasing Classic Aesthetics

NanoPearl Technology – CAD CAM Design makes setup a breeze

Used extensively for All On 4 and Implant Cases

The Pala Mondial tooth line with its wide range of shapes offers a high degree of flexibility in all areas of removable and fixed prosthetics.

Anteriors
RRP \$42.00 ea
Posteriors
RRP \$35.65 ea
Per Card or Less

BUY
30-49
SAVE 10%

BUY 50 +
SAVE 15%

100 + BULK
DISCOUNTS +
INTRO CASE +
CUSTOM BULK
OFFERS
AVAILABLE



Special Intro
Cases
Limited
Available
A2 / A3 or
A3.5 Shades
*New
customer only
– Limit 1

Kulzer Pala PREMIUM 6 / 8 - For The Highest Demands - Mix and Match with Mondial Posteriors

Pala Premium is the tooth line for the highest demands for aesthetics, function, and natural appearance.

Featuring Kulzer nano pearl technology

Removeable and Fixed Implant Retained Overdentures

USED EXTENSIVELY FOR

“IMPLANT RETAINED HIGH-END OVERDENTURES”
PREMIUM ANTERIORS / MONDIAL POSTERIORS

With their fully anatomical morphology and natural light dynamic, Premiums yield the very highest results, particularly for removable and fixed dentures, and thus facilitates effective integration to the residual dentition.

Shades: 16 VitaPan Shades (A1 – D4),
2 Bleach Shades (BL2-BL3)



Anteriors
RRP \$57.25 ea
Per Card or Less

Posteriors
RRP \$47.40ea
Per Card or Less



Valid JULY / AUGUST 2025 - All Prices Exclude GST – Prices Do Not Include Freight - Durodent Reserves The Right To Amend Errors and Prices – Subject to stock availability

ONLINE STORE

www.durodent.com.au

SAVE TIME - ORDER ONLINE



ADD TO CART



Speak to our team

03 9720 6700