

MAY / JUNE 2025 HIGHLIGHTS

durodent dental

PALA DENTURE ACRYLIC RANGE

HEAT CURE

Paladon 65
Strong Heat Cure
Short and Long
Cure Method

SAVE 10%

1kg / 500ml Pack
\$143.28 + GST



SHADE COLOUR MATCH
Heat Cure / Self Cure Shades

SELF CURE / COLD CURE

Pala Dur
Standard Repair - Non
Colour Stable

SAVE 10%

1kg / 500ml Pack
\$221.58 + GST

Pala Press
Colour Stable
All Indications

SAVE 10%

1kg / 500ml Pack
\$268.47 + GST



**R50 Veined
#4**

**Pink Veined
#3**

Clear

*Colour shown on print or
screen may have variations



Pala Xpress
Colour Stable
Longer Working Time
All Indications - Injection

SAVE 10%

1kg / 500ml Pack
\$332.10 + GST

Pala Xtreme
(High Impact Self Cure)
Colour Stable - All Indications

SAVE 10%

1kg / 500ml Pack
\$373.05 + GST

MULTIBUY PACKS
EXTRA DISCOUNTS
6KG + 10KG PACKS

*Prices based on 1kg / 500ml Pack

EOFY

EQUIPMENT

Upgrading Equipment This EOFY Tax Time??

Contact Us For EOFY Equipment Offers

You may be eligible for the small business instant asset write-off up to
\$20,000 - Speak to your Tax Advisor



WASSERMANN
Hamburg - Germany



Renfert
making work easy



OMEC



MOTION



al dente
dentalprodukte



SAESHIN



KULZER
MITSUI CHEMICALS GROUP



**KEYSTONE
INDUSTRIES**

Subscribe @ Durodent.com.au for our monthly email broadcast
Keep informed – New products, highlights, special offers

SUBSCRIBE

Valid MAY / JUNE 2025 - All Prices Exclude GST – Prices Do Not Include Freight - Durodent Reserves The Right To Amend Errors and Prices – Subject to stock availability

ONLINE STORE
www.durodent.com.au

SAVE TIME - ORDER ONLINE



Speak to our team
03 9720 6700

LAB PUTTY

Zhermack Zetalabor Lab Putty - 80 Shore A

C-Silicone Lab Putty
Market Leader
Heat Resistant 140°C
High Elastic Resistance

10% OFF



2.6KG + 1 ACT

\$79.20

ex GST

5KG + 2 ACT

\$121.50

ex GST

LC SPECIAL TRAY

Ernst Hinrich Hinritray LC Tray 50pc Box

Light-curing tray material for
fabricating functional trays,
custom trays, bite-blocks,
registration bite plates and for
other uses in the laboratory.



PINK or BLUE

\$99.90

ex GST

**BUY 3
SAVE 10%**

TRANSPARENT

\$110.00

ex GST

**BUY 6
GET 1 FREE**



GINGIVAL MASK

Zhermack Gingifast Standard or CAD

Elastic 40 Shore A

2 x 50ml

CAD Version is SCANNABLE

includes 10ml Separator and
Mixing Tips

*Stock Expiry Dates 12/25
(Standard) & 07/25 (CAD)

20% OFF



STANDARD*

\$60.00*

ex GST

CAD SCAN*

\$71.20*

ex GST

ALGINATE

Kulzer Alginoplast Fast Set 500g

High Precision Alginate

Our Top Selling Alginate

Industry Favourite



\$25.90

ex GST

**BUY 5 BAGS
5% OFF**

\$24.60ea + GST

**BUY 10 BAGS
10% OFF**

\$23.31ea + GST

**BUY 20 BAGS
20% OFF**

\$20.72ea + GST

ALUMINIUM OXIDE

Ernst Hinrichs Edelkorund Aluminium Oxide

Available in 50um (Fine), 110um (Medium) or 250um (Coarse)

10% OFF



5KG TUB

\$53.10

ex GST



20KG BOX

\$190.80

ex GST

SPECIALTY STONE

Garreco Express Set Stone Type IV Pink - Box 10kg

10% OFF

Exceptionally fast
setting stone with an
early high
compressive
strength, which
allows separation of
the cast from the
impression in 5
minutes.

Dry PSI 14,000
Set Time: 2 Mins
Exp: 0.15%
Ratio: 23ml/100g



\$126.00

ex GST

Valid MAY / JUNE 2025 - All Prices Exclude GST - Prices Do Not Include Freight - Durodent Reserves The Right To Amend Errors and Prices - Subject to stock availability

ONLINE STORE

www.durodent.com.au

SAVE TIME - ORDER ONLINE



ADD TO CART



Speak to our team

03 9720 6700

Hammacher Instrumente

Quality Stainless Steel Lab Instruments
and Pliers Handmade In Germany

10% OFF
ALL IN STOCK
HAMMACHER
ITEMS



PMMA Milling Discs For Premium Denture Base and Teeth

**Dentsply
Sirona**

Multi-Layer PMMA Tooth Shades

Full Shade Range
13 Intricate Layers

Resulting in highly esthetic
outcomes, dentin and enamel
blended seamlessly.



20mm x 98mm

\$111.14 ex GST

BUY 3 + SAVE 10%

Denture Base High Impact PMMA Lucitone

30mm or 25mm x 98mm

Original Pink or Light Pink



30mm x 98mm

\$104.52 ex GST

BUY 3 + SAVE 10%

Laboratory Tungsten Carbide Burs and Cutters



From **\$22.90ea + GST**

DISCOUNT OFFERS

BUY 6 OR MORE

SEE FULL RANGE AND ALL OPTIONS ONLINE AT
WWW.DURODENT.COM.AU

*Can be mixed styles

Standard
Fine
Superfine



Kulzer Signum Opaque F 1 x 3g

Signum opaque F is a light cured single component material used for the colour masking of metal frames.

Available In Shades
A1 – D4
OP (Opaque Pink)
OR (Opaque Red)

*Some shades special order.



10% OFF

\$72.00 ex GST

Vonax Polishing Compound Bar For Acrylic Plastics

Compound for finishing dental plastics.

Polyester Lacquer Finishing Composition

Beige

Approx 730g



10% OFF

\$31.50 ex GST

Kulzer Signum Connector 1 x 5ml 64714211

Bonding Composite To Acrylic

Light curing bonding agent between prosthetic acrylic materials and light curing veneering materials.



10% OFF

\$181.53 ex GST

Kulzer Pala Polish Cream Paste (For Acrylic and Composite) - 80g Jar 66067855

High Shine polishing paste.

Simultaneously removes scratches and leaves a high shine on dental appliances.



10% OFF

\$48.78 ex GST

Polishing Buffs and Brushes

ALL "LATHE" BUFFS AND BRUSHES

10% OFF



Hatho - Renfert - Grobit - Aus Buff - Stitched - Non Stitched - Poly - Slimline - Calico



ADD TO CART

ONLINE ORDER FREIGHT OFFERS

Free Shipping or Subsidised Flat Rate Shipping Offers Available for most orders valued over \$250.00 ex GST with 2 or more items when placed Online through our Web Shop

*Conditions Apply

Excludes some items (eg. Stone, Equipment, Large Items) and delivery locations. It will tell you at checkout if promo applies. Freight promo not available for phone orders.

Kulzer DELARA or MONDIAL Acrylic Teeth

CUSTOM STOCK HOLDING

250 - 348 Cards Per Cabinet

Receive The Following Benefits:

Fully Personalised Shade and Mould Selection

25% DISCOUNT OFF SETUP TEETH

1 x FREE 8 Draw Wooden Labelled Cabinet

2 x FREE Shade Guides

On Going Tooth Discount 10% Off Minimum

***Payment Plan Options** subject to Credit Terms



KULZER
MITSUI CHEMICALS GROUP

DELARA SETUP FROM
\$2806.00 GST FREE

MONDIAL SETUP FROM
\$7175.00 GST FREE



DENTURE TEETH

Kulzer DELARA 6/8
Small Assortment Intro Case
66077520

LIMITED OFFER - NEW USERS

EQUALS 40% DISCOUNT OFF RRP

**ALLOWS YOU TO GIVE DELARA
A GOOD GO AT GREAT PRICE**

\$753.00 ex GST

Includes Shades A2/A3
48 Anteriors / 36 Posteriors
Shade Guide
All Moulds



SPECIAL OFFER

DENTURE TEETH

Kulzer MONDIAL 6/8
Small Assortment Intro Case

LIMITED OFFER - NEW USERS

**EQUALS 40% DISCOUNT OFF RRP
GIVE THEM A GO AT GREAT PRICE**

\$1913.00 ex GST

48 Anteriors / 36 Posteriors - 2 Draws

SELECTION OPTIONS:
A2/A3 - 66041207
A2/A3/A3.5 - 66059081
A3/A3.5 - 66041206



SPECIAL OFFER

ONLINE SHOP CLEARANCE OFFERS

OVERSTOCKED
END OF LINE
SHORT CODE
DISCONTINUED

www.durodent.com.au/clearance



Valid MAY / JUNE 2025 - All Prices Exclude GST - Prices Do Not Include Freight - Durodent Reserves The Right To Amend Errors and Prices - Subject to stock availability

ONLINE STORE
www.durodent.com.au

SAVE TIME - ORDER ONLINE



Speak to our team
03 9720 6700

DURODENT ACRYLIC TEETH RANGE

durodent
dental



KULZER
MITSUBISHI CHEMICALS GROUP

DELARA 3 Layer Modern Aesthetics

Everybody's Delara – a modern all-round tooth line for everyday use with lifelike aesthetics.

A single tooth line that meets every need of dental laboratories – delivering aesthetics, form and function, with a lively surface structure and unique layering, along with modern and bright shading for a fresh and vivid appearance.

3 Layer Incomp Cad Cam

Shades Available: A1 – D4 + 2 Bleach Shades BL2 / BL3

Moulds: 18 Upper / 5 Lower Anterior Moulds - 8 Posteriors Moulds



Only
\$14.85
Per Card or Less

BUY 50-99
SAVE 10%
\$13.36 ea

BUY 100+
SAVE 15%
\$12.62 ea

INTRO CASE +
CUSTOM BULK
OFFERS
AVAILABLE

Kulzer Pala MONDIAL Acrylic Teeth – CAD CAM Designed Perfection – Wear Resistance Tech

Superior Wear Resistance Properties For Fixed or Removable Prosthetics

Highly Pleasing Classic Aesthetics

NanoPearl Technology – CAD CAM Design makes setup a breeze

Used extensively for All On 4 and Implant Cases

The Pala Mondial tooth line with its wide range of shapes offers a high degree of flexibility in all areas of removable and fixed prosthetics.



Special Intro Cases
Limited Available
A2 / A3 or A3.5 Shades
*New customer only –
Limit 1

Anteriors
RRP \$40.60 ea
Posteriors
RRP \$34.45 ea
Per Card or Less

BUY
30-49
SAVE 10%

BUY 50 +
SAVE 15%


100 + BULK
DISCOUNTS +
INTRO CASE +
CUSTOM BULK
OFFERS
AVAILABLE

Kulzer IDEALIS 8 Posteriors – Flat Cusp

5 Deg Flat Cusp Posteriors
Superior Wear Resistance Properties

Ideal posteriors for implants and prostheses for the elderly patients or specific cases requiring low cusps. CAD CAM Manufactured

**Perfectly matched with Kulzer Premium,
Mondial and Delara Anteriors**



Posteriors
RRP \$39.41 ea
Per Card or Less

Kulzer Pala PREMIUM 6 / 8

For The Highest Demands - Mix and Match with Mondial

Pala Premium is the tooth line for the highest demands for aesthetics, function, and natural appearance.

Featuring Kulzer nano pearl technology

Removeable and Fixed Implant Retained Dentures

With their fully anatomical morphology and natural light dynamic, Premiums yield the very highest results, particularly for removable and fixed dentures, and thus facilitates effective integration to the residual dentition.

Shades: 16 VitaPan Shades (A1 – D4), 2 Bleach Shades (BL2-BL3)



Anteriors
RRP \$55.35 ea
Per Card or Less



Posteriors
RRP \$45.83 ea
Per Card or Less

Valid MAY / JUNE 2025 - All Prices Exclude GST – Prices Do Not Include Freight - Durodent Reserves The Right To Amend Errors and Prices – Subject to stock availability

ONLINE STORE

www.durodent.com.au

SAVE TIME - ORDER ONLINE



ADD TO CART



Speak to our team

03 9720 6700