

DuraSplint™

THE NIGHTGUARD WITH FLEX

DuraSplint is the superior and optically clear nightguard with FLEX.

- Optical clarity is nearly invisible while wearing
- Composite formulation, a mixture of powder and liquid, can be press-packed or poured
- Material is stronger with greater stability, ideal for bruxism patients
- Amine-free to never yellow
- MMA-free* to eliminate allergy problems
- Self-adjusting when placed in warm water before inserting in the mouth
- Wearability and patient comfort are rated excellent



Reliance Dental Manufacturing, LLC.

5805 W. 117th Place • Alsip, IL 60803
RelianceDental.net



* MMA- methyl methacrylate



DuraSplint™

Clearly Different

Advantages for both dentist and patient



Dentist



Patient

1. Reliance Dental's flexible **DuraSplint** is self-adjusting. This superb accuracy permits the closest fit.
2. Better fit reduces post-operative chair time.
3. Dimensional stability means fewer, if any, adjustments which are done conventionally.
4. Eliminates allergy problems – no methyl methacrylate like acrylic splints.
5. Ramps and guide planes can be added to reposition the mandible.
6. Chairside repairs can be made with cold cure acrylic.

1. **DuraSplint** is self-adjusting in the mouth when pre-softened in warm water.
2. Better fit reduces post-operative chair time, provides greater comfort, resistance to slipping and rugged wearability.
3. Amine-free – will not yellow with age; always transparent.
4. **DuraSplint** will protect against bruxing and grinding.
5. Breakage-resistant composite resin.

