

# LABO MATE 100 ARTICULATOR



NO.	ITEM	品名	NO.	ITEM	品名
(1)	02092	Main Adjust Screw	(8)	02083	MOUNTING PLATE (WITH MAGNET)
(2)	02093	Main Adjust Nut	(9)	02015	Magnet for Mounting Plate
(3)	02089	Spring	(10)	ART-10022	White iron rod
(4)	02095	Plastic Setting Screw	(11)	ART-10021	STEEL BAR 2PCS/BAG
(5)	02090	Screw for Setting Incisal Pin	(12)	02087-1	Rubber Cover 2pcs/set
(6)	02091	Incisal Pin with Hole	(13)	02085	6x0.6mm steel pan
(7)	02094	Small Guide Pin	(14)		



Metal Mounting Plate  
(without magnet)  
Item NO. 02082  
Plastic Mounting Plate  
(with magnet)  
Item No. 02083



0° 5° 10° 20° 30°  
Incisal Guide Table  
Item No. 02084 complete set  
02084-0°  
02084-5°  
02084-10°  
02084-20°  
02084-30°



02099 02098 02097  
Occlusal Plate Set Include:  
Sphere, plane, stand  
Item No. 02081 complete set  
02099 sphere plate  
02098 plane plate  
02097 stand