

**MADE IN GERMANY**


Wassermann is a medium-sized, owner-operated company with head office in Hamburg. We design and produce high-quality technical devices in Germany. Our main users are dental laboratories, dental clinics and various industrial areas worldwide.

**Why a Wassermann product?**

Wassermann products are manufactured from high-quality materials and the proportion of plastics used is intentionally kept as low as possible. Our devices are known for their suitability for daily use, functionality and longevity. As a result, the users place high levels of trust in the Wassermann brand. We feel it is our duty to support our customers on a daily basis when it comes to safety, service and productivity. Technological perfection, optimum quality and especially sustainability are the benchmarks we set. Thanks to the combination of these values, the family-owned company, since its establishment in 1927, has gained worldwide recognition as a manufacturer of high-quality products.


 WP-Ex 2000 II (2 settings, 1500 / 3000 min<sup>-1</sup>), basic equipment


WP-Ex 2000 II with special accessories daylight LED spot and suction protective grids

**WP-Ex 3000 II and WP-Ex 2000 II**

**Polishing Compact Units**


Further information are available on our website at [www.wassermann.hamburg](http://www.wassermann.hamburg)

 Your Dental Dealer
 

Wassermann Dental-Maschinen GmbH  
 Rudorffweg 15-17 · 21031 Hamburg · Germany  
 Phone: +49 (0)40 730 926-0 · Fax: +49 (0)40 730 37 24  
[info@wassermann-dental.com](mailto:info@wassermann-dental.com) · [www.wassermann.hamburg](http://www.wassermann.hamburg)

# Polishing Compact Units

## WP-Ex 3000 II & WP-Ex 2000 II

Powerful and durable: The WP-Ex 3000 II and WP-Ex 2000 II polishing compact units for daily use at labs guarantee first-class polishing and extraction. Due to their compact size they can be used in the tiniest space. They stand out thanks to their powerful and low-noise extraction and the strong, durable and maintenance-free induction motor. Convenient: The WP-Ex 3000 II is easy to operate due the simple push of a button and the WP-Ex 2000 II even with the ball of the thumb.

WP-Ex 3000 II comes as standard with a clear multifunctional display, a motor which is variable in steps and an adjustable daylight LED spot

Emergency off-switch and restart protection guarantee safety at work

Impact-resistant, adjustable safety screen

Easy to clean housing, made almost entirely of stainless steel, integrated drawer

Numerous storage options

Integrated day-light LEDs

Two easily removable full rubber polishing troughs



### Major Features

- Motor with jacket-cooling, dust protected
- Two micro fine dust filterbags with large capacity which are easy to replace guarantee best filter results
- **Included in delivery:**  
 Brush spindle right, 92 mm / Brush spindle left, 92 mm  
 Micro fine dust filterbag, 1 set = 2 pieces  
 Storage bowl (W 257 x H 40 x D 198 mm)  
 Storage bowl (W 190 x H 17 x D 150 mm)  
 Drawer (W 621 x H 23 x D 219 mm)  
 Polishing trough, 2 pieces  
 Supply line

Accessories	Item no.
Daylight LED spot (WP-Ex 3000 II comes as standard with daylight LED spot)	146390
Suction protective grid, 1 piece	146396
Brush spindle right, 92 mm	152022
Brush spindle left, 92 mm	152025
Brush spindle right, 125 mm	152043
Brush spindle left, 125 mm	152044
Stone holder right	152056
Stone holder left	152057
Drill chuck right, 1-6 mm clamping	152049
Storage bowl (W 257 x H 40 x D 198 mm)	830409
Storage bowl (W 190 x H 17 x D 150 mm)	395004
Spare micro fine dust filterbag, 1 set =2 pieces	611005
<b>WP-Ex 3000 II</b>	<b>Item no. 146979</b>
<b>WP-Ex 2000 II</b>	<b>Item no. 146978</b>
Voltage	220–240 V / 50/60 Hz
Power consumption	2.8 A
Output	650 W
Fan capacity	445 m <sup>3</sup> /h (free-blowing)
W x H x D	680 x 486 x 527 mm (Height with LED spot: + 140 mm)
Weight with motor	52.0 kg
Sound level	≤ 65 dB (A)
Speed	1500 / 3000 min <sup>-1</sup>
Speed WP-Ex 3000 II	additionally adjustable in steps 750–6000 min <sup>-1</sup>