

Delara®

Mould chart

NEW



Delara®

A modern all-round tooth line for all indications with lifelike aesthetics.

Giving a hand to oral health.



KULZER
MITSUI CHEMICALS GROUP

Anteriors

SMALL

T 42 ▽		42,4 9,6 8,0	S 44 □		44,5 10,4 8,3
O 43 ○		44,2 10,5 8,3	T 44 ▽		45,1 10,5 8,2
S 43 □		43,9 9,6 8,2	R 45 □		45,0 10,6 8,2

MEDIUM

T 45 ▽		45,7 11,3 8,4	O 47 ○		48,1 11,4 9,1
T 46 ▽		46,5 10,7 8,7	R 47 □		47,5 11,5 8,8
R 46 □		46,5 10,1 8,8	S 48 □		49,0 11,3 9,1

LARGE

R 48 □		48,8 13,0 9,3	S 50 □		51,0 11,8 9,6
T 49 ▽		49,6 12,2 9,1	O 51 ○		52,2 12,4 9,8
R 49 □		48,7 10,9 9,3	R 52 □		52,1 13,6 9,8

L 33



33,6
9,4
5,2

L 38



36,6
10,5
5,8

L 35



34,6
9,5
5,4

L 39



38,3
10,8
5,9

L 36



35,6
9,1
5,5

Posteriors

30

U

L



29,5
9,1
31,4
8,1

32

U

L



31,0
9,5
33,1
8,8

34

U

L



33,0
10,2
35,3
9,3

36

U

L



34,8
10,9
37,3
9,9

Mould numbers and possible combinations



Upper anteriors
T 42
T 44
T 45
T 46
T 49

Lower anteriors
L 33
L 33 + L 35
L 33 + L 35 + L 36
L 35 + L 36 + L 38
L 38 + L 39

Posteriors
30 + 32
30 + 32
32 + 34
32 + 34
36



O 43
O 47
O 51

L 33 + L 35 + L 36
L 38 + L 39
L 38 + L 39

30 + 32
34
36



R 45
R 46
R 47
R 48
R 49
R 52

L 33 + L 35
L 33 + L 35 + L 36
L 38 + L 39
L 36 + L 38
L 38 + L 39
L 38 + L 39

32
32 + 34
34
34
34 + 36
36




S 43
S 44
S 48
S 50

L 33 + L 35 + L 36
L 33 + L 35 + L 36 + L 38
L 38 + L 39
L 38 + L 39

30 + 32
30 + 32 + 34
34 + 36
36


Examples




T 42
42,4
9,6
8,0



30 U
29,1



L 33
33,6
9,4
5,2



30 L
31,4
8,1

Everybody's Delara



Delara meets every need offering a lively surface structure and unique layering, along with modern and bright shading for a fresh and vivid appearance.

- **Easy handling for predictable results:** freedom during set-up without excessive grinding or changing the basic shape and layering effects. The polishing buffer, basal grinding reserve and optimized incisal edge simplify the set-up process.
- **Flexibility:** bigger interdental space, full range of shades, process-optimized and anatomical forms – for all indications and set-up options.
- **Aesthetics:** blueish effect for a fresh and lively appearance, along with modern shading, a vivid surface structure and dynamic light refraction.

That's why we call it
"Everybody's Delara".

“Delara anterior teeth look so natural and beautiful in the mouth thanks to the unique surface structure. In addition, they are easy to set up, especially the posterior teeth. With Delara you have the freedom of set-up with almost no grinding necessary.”



Ilkka Garaisi

As dentist (clinical dental technician) he has over 25 years of experience and treated over 7.000 patients with full dentures. He is president of EHT-federation in Finland, teacher of denturists, and representative of Finland in the IFD.





Manufacturer

Kulzer GmbH
Schlätterstraße 2
88142 Wasserburg, Germany

Contact

Kulzer Australia Pty Ltd
info.australia@kulzer-dental.com

kulzer.com/delara